


# Oriel® VeraSol™ Solar Simulator\*



- Output beam size: 2 in. x 2 in. (51 mm x 51 mm)
- Factory certified IEC, ASTM, JIS AAA rated; CE certified
- Variable output adjustment from 0.1 to 1.1 SUN
- Fast turn on time, <100 ms - no shutter required
- 10,000 hour LED lifetime - no bulb replacement required
- User settable spectral control
- User adjustable calibration
- Flexible mounting orientation
- PV cell placement indicator
- USB LED control

\* Patent Pending

Newport's Oriel® brand is proud to introduce the first in a series of innovative solar simulator designs utilizing LED light sources. The system is made up of the LSS-7120 Controller and the revolutionary LSH-7520 LED Source Head. The system provides a variable output from 0.1 SUN to 1.1 SUN over a 2 inch x 2 inch (51 mm x 51 mm) illumination area.

The LSH-7520 is certified AAA rated to IEC 60904-9, JIS C8904-9, and ASTM-E927-10 for Spectral Match, Non-Uniformity of Irradiance, and Temporal Stability. Oriel puts each LED based system through rigorous testing for all aspects of the standards to ensure compliance and provides a certificate of calibration.

This innovative design independently drives multiple LEDs at 19 individual wavelengths spaced over the spectrum from 400 nm to 1100 nm to ensure a spectral match meeting Class A certification. The independent control of each LED allows the user to tune the output to the specific requirements of the test. This adjustment can be done at the front panel for the entire spectrum or individual band adjustment. The user can adjust individual LED output through the USB port. These spectral matches can be saved as a preset value in the control instrument.

The typical lifetime of LEDs of over 10,000 hours removes the need to change expensive and hazardous bulbs. Additionally, the LEDs come up to a stable operating power within 100 msec. This allows the unit to be shut off between tests without having to wait for the lamp to stabilize or sacrificing the integrity of the test results. This results in an effective lamp usage time several times over a conventional lamp and eliminates the need for shutters on the output.

The use of solid state LEDs and robust optical mounting hardware makes it possible to orient the head in any position to fit your application needs. To aid in positioning the light source, a variable height adjustment mount is included along with a laser pointer indicator to ensure your test cell is always at the correct placement under the head. The flexibility of the VeraSol™ Solar Simulator allows the user to tailor the simulator to the test and not the test to the solar simulator.


Solar simulation has entered a new era. We believe you will continue to find features and benefits of the VeraSol™ that will improve your testing capabilities.

**Oriel**  
**INSTRUMENTS**  
A Newport Company

# Oriel® VeraSol™ Solar Simulator\*




## LSH-7520 LIGHT SOURCE SPECIFICATIONS

Illumination Area Maximum	2 inches x 2 inches (50 mm x 50 mm)
Calibrated Output Power	100 mW/cm <sup>2</sup> (1.0 SUN)
Wavelength Range	400 nm - 1100 nm
Uniformity Classification <sup>1</sup>	A - IEC 60904-9 2007, JIS C8904-9 2017, ASTM E927-10 2015
Temporal Stability	A - IEC 60904-9 2007, JIS C8904-9 2017, ASTM E927-10 2015
Spectral Match	A - IEC 60904-9 2007, JIS C8904-9 2017, ASTM E927-10 2015
Source Orientation <sup>2</sup>	0°, 180°
Weight (head only)	17.7 kg (39.0 lbs)
Nominal Working Distance	8.0" (200 mm)
Height Adjustment	8.9" - 14.1" (226 mm - 358 mm)



## LSS-7120 LED SOLAR SIMULATOR SPECIFICATIONS

Independent Band Control	6 Bands (400 - 500 nm, 500 - 600 nm, 600 - 700 nm, 700 - 800 nm, 800 - 900 nm, 900 - 1100 nm)
Variable Output Control <sup>3</sup>	0.1 to 1.1 SUN
Auxiliary Functions	10 user settable presets User settable output calibration LED fault detection
I/O Connectors	USB 2.0 (B-type) DB-60; LSH interconnect port
Power Requirements	100-240 VAC, 50-60 Hz, 300W
Size	4.0" x 8.5" x 14.0" (102 mm x 216 mm x 356 mm)
Weight	4.9 kg (10.8 lbs)
Operating Temp Range	20°C to 30°C
Storage Temp Range	-40°C to 70°C
Humidity	<85%, relative, non-condensing
Compliance	CE



Typical Results



**Newport Corporation**  
31950 Frontage Road Bozeman, MT, USA

PHONE: 1-800-459-9459 1-406-556-2481 FAX: 1-406-586-9405 EMAIL: [OrielPV.sales@newport.com](mailto:OrielPV.sales@newport.com)  
Complete listings for all global office locations are available online at [www.newport.com/contact](http://www.newport.com/contact)

[www.newport.com](http://www.newport.com)

	PHONE
Belgium	+32-(0)0800-11 257
China	+86-10-6267-0065
France	+33-(0)1-60-91-68-68
Japan	+81-3-3794-5511
Taiwan	+886-(0)2-2508-4977

	EMAIL
Belgium	<a href="mailto:belgium@newport.com">belgium@newport.com</a>
China	<a href="mailto:china@newport.com">china@newport.com</a>
France	<a href="mailto:france@newport.com">france@newport.com</a>
Japan	<a href="mailto:spectra-physics@splasers.co.jp">spectra-physics@splasers.co.jp</a>
Taiwan	<a href="mailto:sales@newport.com.tw">sales@newport.com.tw</a>

	PHONE
Irvine, CA, USA	+1-800-222-6440
Netherlands	+31-(0)30 6592111
United Kingdom	+44-1235-432-710
Germany / Austria / Switzerland	+49-(0)6151-708-0

	EMAIL
Irvine, CA, USA	<a href="mailto:sales@newport.com">sales@newport.com</a>
Netherlands	<a href="mailto:netherlands@newport.com">netherlands@newport.com</a>
United Kingdom	<a href="mailto:uk@newport.com">uk@newport.com</a>
Germany / Austria / Switzerland	<a href="mailto:germany@newport.com">germany@newport.com</a>

Newport Corporation, Irvine and Santa Clara, California and Franklin, Massachusetts; Evry and Beaune-La-Rolande, France; Stahnsdorf, Germany and Wuxi, China have all been certified compliant with ISO 9001 by the British Standards Institution.

1. Certified over a 1.9 inch (48 mm) diameter circle within a 2 x 2 inch illumination area.
2. A mount plate can be added to allow for 90° orientations.
3. Output intensity (other than 1.0 SUN, which is calibrated at the factory) is approximate and Class AAA performance is only certified at 1.0 SUN.

VERASOL REV06 082919