

USER'S GUIDE

DC-250MHz Electro-Optic Phase Modulators

Models 400X, 406X

U.S. Patent #5,189,547



3635 Peterson Way • Santa Clara, CA 95054 • USA
phone: (408) 980-5903 • fax: (408) 987-3178
e-mail: techsupport@newfocus.com • www.newfocus.com

Warranty

Newport Corporation guarantees its products to be free of defects for one year from the date of shipment. This is in lieu of all other guarantees, expressed or implied, and does not cover incidental or consequential loss.

Information in this document is subject to change without notice.

Copyright 2012, 2001-1998, Newport Corporation. All rights reserved.

The New Focus logo and symbol are registered trademarks of Newport Corporation

Document Number 400412 Rev. J

Contents

Operation

Introduction	5
Quick Start	7
Using the Modulator	8

Principles of Operation

The Electro-Optic Effect	13
Creating Sidebands	15

Characteristics

Specification	16
Definitions of Specifications	18

Customer Service

Technical Support	21
Service	21

Performance Data 22

Operation

Introduction

The New Focus model 400X and 406X series of electro-optic phase modulators provides an efficient means for optical phase modulation in the DC to 250 MHz frequency range. These versatile phase modulators can be operated with low to moderate drive voltages, and their 2-mm apertures make them compatible with most laser sources. Other features of these devices include their high modulation frequency, low insertion loss, good RF shielding, and high optical-power handling capability. The model 400X and 406X series consists of modulators classified according to operating wavelengths, optical-damage thresholds and electronic drive requirements. The operating wavelengths are determined by the anti-reflection coating applied to the surfaces of the electro-optic crystals. Three standard wavelength ranges are offered: 360-500nm, 500-900nm and 900-1600nm. Regarding electronic drive requirements two types of modulators are offered; broadband and resonant.

- Models 4002, 4004 and 4006 Broadband Phase Modulators utilize MgO-doped lithium niobate crystals and can be operated at any frequency from DC to 100 MHz.
- Models 4062 and 4064 KTP Broadband Phase Modulators can be operated at any frequency from 0.01 MHz to 250 MHz.

- Models 4001, 4003 and 4005 Resonant Phase Modulators utilize MgO-doped lithium niobate crystals and operate at a single user-specified frequency anywhere in the range 0.01 to 250 MHz. These devices can only be operated at their resonant frequency but require much lower drive voltages.
- Models 4061 and 4063 Resonant Phase Modulators utilize KTP crystals and operate at a single user-specified frequency anywhere in the range of 0.1 to 250 MHz. These devices can only be operated at their resonant frequency but require much lower drive voltages.

The broadband phase modulators are appropriate for applications where modulation over a broad frequency range is required. For applications requiring phase modulation at a single frequency, the resonant phase modulators are preferred because much higher phase shifts can be achieved with a given drive voltage.


All of these devices employ electro-optic crystals where an applied electric field induces a change in the crystal's refractive index. The Model 400X series phase modulators use magnesium-oxide-doped lithium niobate (MgO:LiNbO₃) crystals as the electro-optic medium. MgO:LiNbO₃ is nonhygroscopic, and so these modulators can be left on an optical table for indefinite periods without requiring a sealed enclosure.

Model 406X series phase modulators use Potassium Titanyl Phosphate (KTP) as the electro-optic medium. Compared with lithium niobate, KTP crystals have higher optical-damage thresholds and reduced optical wavefront distortion. However, KTP crystals require slightly higher drive voltages.

Performance specifications for the 400X series phase modulators are listed in Table 1 on page 16.

Performance specifications for the 406X series phase modulator are listed in Table 2 on page 17.

**Figure 1:
Mechanical
views of the
Model 400X
and 406X**



Quick Start

This section presents a brief introduction to using your high-frequency phase modulator.

- Align a collimated optical beam through the mechanical apertures of the modulator. For all Models 400X and 406X the beam should be polarized vertically with respect to the modulator casing.



Be careful not to exceed the maximum recommended optical power, or damage to the electro-optic crystal could result. (See page 11 for a discussion of optical damage.)

- Drive the modulator with a 50- Ω RF driver tuned to the modulator's resonant frequency. RF powers from 0.1 to 0.5 watts should be sufficient to allow observation of sidebands. Generally, an optical spectrum analyzer with suitable finesse and free spectral range is used to observe the modulation sidebands.



To prevent damaging the electro-optic crystal, do not exceed the modulator's maximum RF drive power (see Table 1 on page 16 and Table 2 on page 17 for maximum RF power).



If the modulator is not driven at or close to its resonant frequency, most of the RF drive power will be reflected, which could cause damage to the driver.

Using the Modulator

When used properly, the 400X and 406X series phase modulators can provide efficient phase modulation with extremely low unwanted amplitude modulation and insertion loss. The key to obtaining this pure phase modulation is good optical alignment of the beam to the crystal's propagation axis, and accurate polarization orientation of the laser's electric field with the crystal's electro-optically active axis.

Aligning an Optical Beam Through the Modulator

To align the module to the optical beam:


1. Use the 1/4-20 (M6 for metric versions) tapped hole located on the base of the module to mount it on an adjustment-positioning device for alignment. We recommend the New Focus Model 9071, 9071M 9081, or 9081M tilt aligner because of its tilt and translation capabilities.
2. Turn on the optical beam, and orient the beam so it is vertically polarized on the input aperture (see Figure 2). It is important to carefully align the

polarization since the crystals used by New Focus are cut so that the beam propagates along the y-axis of the crystal. This orientation minimizes the effects of acoustic resonances but makes it critical that the optical beam be linearly polarized vertically along the z-axis. If the polarization of the optical beam is not properly aligned, the modulators will impose a polarization rotation as well as a phase modulation which can lead to unwanted amplitude modulation if the modulator is followed by any polarizing optics.

3. Position and align the module so that the beam passes through the 2-mm input and output apertures, clearing them without clipping. The beam to be modulated should be collimated with a waist size of less than 2 mm, and such that the Rayleigh range is at least the length of the crystal. To avoid damage resulting from excessive optical intensity, the optical power should be kept below the damage thresholds listed in Table 1. Typically, a good beam size is 250–500 μm .

Figure 2 shows a phase modulator being driven by a source tuned to f_R . The module is mounted on a Model 9071 tilt aligner, and the input beam is vertically aligned.

Figure 2:
Model 400X
and 406X
with vertically
aligned input
beam



Setting Up the Electrical Modulating Input Signal

Using an SMA cable, connect the SMA jack on the modulator to a modulating source appropriate for the type of modulator (resonant or broadband) you are using.

Note

Since the optical alignment of the modulator can be disturbed by the SMA cable, ensure that the SMA orientation is not obstructing the alignment and use a strain relief on the cable.

The Models 4001, 4003, 4005, 4061 and 4063 resonant modulators are tuned to a specific frequency and require very low drive voltages, such as that from a simple crystal oscillator or a function generator that has an output impedance near 50Ω . The reflected power from resonant phase modulators at their specified resonant frequency is very low compared to broadband modulators.

The Models 4002, 4004, 4006, 4062 and 4064, broadband modulators require large drive voltages and have a bandwidth dependent on the impedance of the modulating source. The source must be able to drive an open circuit without causing damage to the source.

Preventing Photorefractive Damage

The electro-optic crystals used in these modulators are susceptible to optical damage through the photorefractive effect. Photorefractive damage becomes more of a problem for high optical powers, short wavelengths, and tightly focused beams. Table 1 on page 16 and Table 2 on page 17 list maximum average power levels above which optical damage will occur.

The photorefractive damage process can occur gradually over days or hours, or, for high optical powers and short wavelengths, this effect can occur over seconds. A damaged crystal will distort a beam, usually by elongating it in the vertical direction. If operating close to the damage threshold, it is a good idea to monitor the transmitted beam periodically for indications of optical damage. Photorefractive damage can be partially reversed by carefully annealing the crystal; please contact New Focus for details on this procedure.

Principles of Operation

The Electro-Optic Effect

Operation of New Focus electro-optic phase modulators is based on the linear electro-optic (or Pockels) effect—the linear dependence of the refractive index on the applied electric field. The effect of an applied electric field on a crystal's refractive index is described by a third-rank tensor r_{ij} . The induced refractive index change caused by an external electric field has the form

$$\Delta n = \frac{1}{2} n_e^3 r_{33} E$$

where Δn is the change in the index of refraction, n_e is the unperturbed index of refraction, r_{33} is the appropriate element in the electro-optic tensor, and E is the applied electric field. We use magnesium-oxide-doped lithium niobate (MgO:LiNbO₃) and Potassium Titanyl Phosphate (KTP) in these modulators because they have a wide spectral-transparency window, a large coefficient, and low RF losses.

These phase modulators consist of an electro-optic crystal of length l with electrodes separated by the crystal thickness d . The electric field is applied along a crystal axis transverse to the direction of optical propagation. Modulation is induced onto the laser beam by aligning the polarization of the input beam with the crystal axis along which the electric field is applied. An electronic signal is then directly modulated onto the laser beam through the electro-optic effect.

For the Models 4002, 4004, 4006, 4062, and 4064 broadband phase modulators, the input electronic signal is applied directly across the crystal's electrodes. So, the optical phase shift obtained by applying a voltage V at the input SMA connector is

$$\Delta\phi = \frac{2\pi}{\lambda} \left(\frac{1}{2} n_e^3 r_{33} \right) \frac{1}{d} V$$

where λ is the free-space wavelength. A commonly used figure of merit for electro-optic modulators is the half-wave voltage, V_π , which is defined as the voltage required to produce a phase shift of 180° . Substituting into the preceding equation yields

$$V_\pi = \frac{\lambda d}{n_e^3 r_{33}}$$

For the Lithium Niobate-based broadband phase modulators, such as Models 4002, 4004 and 4006, V_π is typically 210 volts at $1.06 \mu\text{m}$, corresponding to a modulation depth, β , of 0.015 radians/volt. Note that at other wavelengths these values change proportionately. So, at 532 nm V_π is 105 volts, and β is 0.03 radians/volt.

For KTP-based broadband phase modulators, such as Models 4062, 4064 and 4066, V_π is slightly higher, typically 230 Volts at $1.06 \mu\text{m}$ with proportional dependence on wavelength similar to Models 4061 and 4063.

For Models 4001, 4003, 4005, 4061 and 4063, resonant phase modulators, the crystal is combined with an inductor to form a resonant tank circuit. On resonance, the circuit looks like a resistor whose value depends on the inductor's losses. A transformer is used to match this resistance to the $50\text{-}\Omega$ driving impedance. Putting the crystal in this resonant circuit results in a voltage across the crystal electrodes that can be more than ten times the input voltage across the SMA connector. This leads to reduced half-wave voltages and larger modulation depths when compared to the broadband modulators.

For the resonant phase modulators, the peak phase shift obtained by applying a sinusoidal signal of average power P at the input SMA connector is

$$\Delta\phi = \frac{2\pi}{\lambda} \left(\frac{1}{2} n_e^3 r_{33} \right) \frac{1}{d} \sqrt{2PQ} \sqrt{\frac{L}{C}}$$

where Q is the quality factor of the tank circuit, L is the inductance, and C is the crystal capacitance. For the resonant phase modulators $V\pi$ is typically 16 volts at $1.06 \mu\text{m}$, corresponding to a modulation depth β of 0.2 radians/volt.

Creating Sidebands

A sinusoidal waveform applied to the modulator will produce frequency sidebands which are separated from the optical carrier by the modulation frequency. These modulation sidebands can be observed by looking at the beam with an optical spectrum analyzer. Given an induced peak optical phase shift of $\Delta\phi$ (in radians), the fraction of power transferred to each of the first-order sidebands is $[J_1(\Delta\phi)]^2$, where J_1 is the Bessel function of order 1. The fraction of power that remains in the carrier is $[J_0(\Delta\phi)]^2$, where J_0 is the Bessel function of order 0.

For example, imposing a sinusoidal phase shift of peak amplitude 1 radian to a cw laser beam will transfer 19% of the initial carrier power to each of the first-order sidebands and leave 59% of the power in the carrier. Note that the maximum power that can be transferred to the first-order sidebands is about 34%, and this requires a peak phase shift of 1.8 radians. For the broadband phase modulator, at $1.06 \mu\text{m}$ a 1.8 radian phase shift is achieved with about 120 volts. For the resonant phase modulator, a 1.8 radian phase shift requires an input peak voltage of about 9 volts (0.81 W average power).

Characteristics

Specifications

Table 1: Model 400X Series Specifications

	4001 4003 4005	4002 4004 4006
Wavelength	360-500 nm (4005) 500-900 nm (4001) 900-1600 nm (4003)	360-500 nm (4006) 500-900 nm (4002) 900-1600 nm (4004)
Type	Resonant PM	Broadband PM
Operating Frequency*	.01 - 250 MHz (single frequency)	DC-100MHz
Modulation Depth	0.27 - 0.80 rad/V (364 nm) 0.2 - 0.6 rad/V (532 nm) 0.1 - 0.3 rad/V (1000 nm)	40 mrad/V (364 nm) 30 mrad/V (532 nm) 15 mrad/V (1000 nm)
Max $V\pi$	3.8 - 11.7 V (364 nm) 5 - 16 V (532 nm) 10 - 31 V (1000 nm)	79 V (364 nm) 105 V (532 nm) 210 V (1000 nm)
Material	MgO:LiNbO ₃	MgO:LiNbO ₃
Max Optical Intensity**	0.1 W/mm ² (364 nm) 2 W/mm ² (532 nm) 4 W/mm ² (1064 nm)	0.1 W/mm ² (364 nm) 2 W/mm ² (532 nm) 4 W/mm ² (1064 nm)
Aperture	2 mm x 2 mm	2 mm x 2 mm
RF Bandwidth	2-4 % of center frequency	100 MHz
Connector	SMA	SMA
Impedance	50 Ω	20 pF
Max. RF Power	1 W	10 W
Max VSWR	1.5	N/A

*Resonant modulators must be specified to a single frequency.

**In a 1 mm beam.

Table 2: Model 406X Series Specifications

	4061 4063	4062 4064
Wavelength	500-900 nm (4061) 1000-1600 nm (4063)	500-900 nm (4062) 1000-1600 nm (4064)
Type	Resonant PM	Broadband PM
Operating Frequency*	.01 - 250 MHz (single frequency)	DC - 250 MHz
Modulation Depth	0.16 - 0.3 rad/V (532 nm) 0.08 - 0.16 rad/V (1000 nm)	26 mrad/V (532 nm) 13 mrad/V (1000 nm)
Max $V\pi$	10 - 31 V (532 nm) 19 - 58 V (1000 nm)	115 V (532 nm) 230 V (1000 nm)
Material	KTP	KTP
Max Optical Intensity**	10 W/mm ² (532 nm) 20 W/mm ² (1064 nm)	10 W/mm ² (532 nm) 20 W/mm ² (1064 nm)
Aperture	2 mm x 2 mm	2 mm x 2 mm
RF Bandwidth	2-4 % of center frequency	250 MHz
Connector	SMA	SMA
Impedance	50 Ω	10 pF
Max. RF Power	2 W	10 W
Max VSWR	1.5	N/A

*Resonant modulators must be specified to a single frequency.

**In a 1 mm beam.

Definitions of Specifications

Optical Throughput

Optical throughput is determined by the absorption and scatter in the electro-optic crystal, and by the quality of the anti-reflection coatings on the end faces. Low optical losses are critical in applications of the New Focus phase modulators, so great care is taken to ensure insertion loss is minimized.

Modulation Depth

This describes the magnitude of the phase modulation imposed on the input laser beam by the modulator. This depth is optimized by alignment of the input beam's polarization with the crystal active axis.

Residual Amplitude Modulation (RAM)

RAM is noise in a phase modulation system, and therefore must be minimized. Quality of crystal growth and excellence in the finished crystal design and manufacturing are essential to the elimination of RAM. Experimentally, precise alignment is required to prevent RAM from occurring.

Return Loss

The return loss indicates the quality of impedance matching between the driving source and a resonant phase modulator. Resonant modulators are designed to have impedances very close to 50Ω at resonance, and a high return loss indicates a good impedance match between the driving source and the modulator. With a high return loss, power transfer to the modulator is optimized, and reflected power, which can harm RF drivers, is minimized.

The return loss indicates the fraction of RF power reflected from the modulator back to the driver when the modulator is driven at its resonant frequency. For

a power reflection coefficient R, the return loss in dB is $-10 \log R$. All New Focus resonant phase modulators are tested by measuring return loss versus frequency around the modulation frequency. The results of this test are provided in the specifications table above.

Voltage Standing Wave Ratio (VSWR)

VSWR, like return loss, is another way to specify the quality of impedance matching between RF driver and resonant modulator. It is defined as the voltage ratio between the maximum and minimum of the standing wave that occurs because of an impedance mismatch. A VSWR value of 1 indicates a perfectly matched system. Given a return loss RL (in dB), the VSWR can be found from

$$\text{VSWR} = \frac{1+10^{(-RL)/20}}{1-10^{(-RL)/20}}$$

Customer Service

Technical Support

Information and advice about the operation of any New Focus product is available from our applications engineers. For quickest response, ask for “Technical Support” and know the model number and serial number for your product.

Hours: 8:00–5:00 PST, Monday through Friday (excluding holidays).

Toll Free: 1-866-NUFOCUS (1-866-683-6287) (from the USA & Canada only)

Phone: (408) 980-5903

Support is also available by fax and email:

Fax: (408) 987-3178

Email: techsupport@newfocus.com

We typically respond to faxes and email within one business day.

Service

In the event that your modulator malfunctions or becomes damaged, please contact New Focus for a return authorization number and instructions on shipping the unit back for evaluation and repair.

Performance Data

Model Number: _____

Serial Number: _____

Frequency: _____

Wavelength: _____

Input RF Power: _____

Return Loss: _____

VSWR: _____

Q: _____

