

Broadband Ultrafast Variable Attenuators (VA-BB)



Motorized VA-BB Broadband Ultrafast Variable Attenuator.
(Post and post holder sold separately)



Manual VA-BB Broadband Ultrafast Variable Attenuator

Features

- Compact, robust design comes in manual or motorized version
- Wavelength ranges of 400 – 700 nm, 690 – 1040 nm or 1000 – 1600 nm
- High contrast attenuation (1000:1)
- Minimal stretching for ultrafast pulses >70 fs
- High optical damage threshold over 10 mm clear aperture
- P-polarized collinear output

Newport's broadband ultrafast variable attenuators (VA-BB) provide continuous attenuation over a series of broad wavelength ranges. The VA-BB series is designed to provide high extinction ratio attenuation of linearly polarized fundamental and harmonic frequencies of light for most ultrafast oscillators. The VA-BB series is available in manual or motorized versions. The motorized version can be used in open or closed-loop operation.

The VA-BB series is based on zero-order half-waveplates and Glan-laser polarizers. For more detail on the principle of operation of these devices, please refer to [Newport's Application Note 26 – Variable Attenuator for Lasers](#). These attenuators are designed to handle average powers up to 100 Watts (see table on opposite side of this datasheet). Specified optics sets are optimized for wavelength ranges of 400 – 700 nm, 690 – 1040 nm or 1000 – 1600 nm.

The optics and high-precision opto-mechanics are incorporated into a compact housing, flexibly designed for post mounting at varying optical axis heights. P-polarized light is transmitted collinear to the input beam path with minimal loss while S-polarized light is rejected into a beam dump on the side of the housing. M-4 and 8-32 mounting threads are included on the bottom and one side of the housing. This allows for horizontally and vertically polarized output.

The manual versions of the VA-BB series come with a manual rotary stage containing the half-waveplate. The user manually rotates the stage to control attenuation.

The motorized versions of the VA-BB series come with a motorized rotary stage containing the half-waveplate. Software is included which allows the user to calibrate the VA-BB at a single wavelength. When using Newport family power meters and detectors (purchased separately), this calibration is automated. When using third-party power meters and detectors, software prompts will step the user through building up the calibration curve. The calibration curve is created by rotating the half-waveplate incrementally through a full attenuation cycle. The absolute output power at each increment is measured with the power meter and detector, and a calibration curve is fit to the data points.


Once the calibration curve is created, the motorized VA-BB can be used in open-loop mode. The power meter and detector are no longer needed as the desired output power is specified as an absolute or percentage value. The half-waveplate is automatically rotated to the appropriate position matching that value on the calibration curve.

The motorized VA-BB can also be used in closed-loop mode if desired. This option requires a beam sampler assembly (purchased separately) to pick off a portion of the transmitted beam for real-time monitoring by a Newport power meter and detector. The closed-loop control will fine adjust to the position of the half-waveplate to reach the requested absolute power.




Broadband Ultrafast Variable Attenuators (VA-BB)

Open Loop Configuration



Closed Loop Configuration




Maximum laser power vs. beam diameter

Beam diameter (mm)	5	4	3	2	1
Laser Power (W)	100	64	36	16	4

Typical Specifications

P-polarization attenuation	1000:1
Clear aperture	10 mm diameter
Damage threshold	500 W/cm ² CW, 2 J/cm ² with 8 nsec pulses at 1064 nm, typical



Newport Corporation, Global Headquarters
1791 Deere Avenue, Irvine, CA 92606, USA

www.newport.com

PHONE: 1-800-222-6440 1-949-863-3144 FAX: 1-949-253-1680 EMAIL: sales@newport.com
Complete listings for all global office locations are available online at www.newport.com/contact

PHONE

Belgium +32-(0)800-11 257
China +86-10-6267-0065
France +33-(0)1-60-91-68-68
Japan +81-3-3794-5511
Taiwan +886 -(0)2-2508-4977

EMAIL

belgium@newport.com
china@newport.com
france@newport.com
spectra-physics@splasers.co.jp
sales@newport.com.tw

PHONE

Irvine, CA, USA +1-800-222-6440
Netherlands +31-(0)30 6592111
United Kingdom +44-1235-432-710
Germany / Austria / Switzerland +49-(0)6151-708-0

EMAIL

sales@newport.com
netherlands@newport.com
uk@newport.com
germany@newport.com

Newport Corporation, Irvine and Santa Clara, California and Franklin, Massachusetts;
Evry and Beaune-La-Rolande, France; Stahnsdorf, Germany and Wuxi, China have
all been certified compliant with ISO 9001 by the British Standards Institution.

DS-101203