

# CONEX-AGAP

## Agilis-D Controller with Strain Gages Feedback

---



 **Newport**<sup>®</sup>

**Command Interface  
Manual**

V2.0.x

©2019 by Newport Corporation, Irvine, CA. All rights reserved.

Original instructions.

No part of this document may be reproduced or copied without the prior written approval of Newport Corporation. This document is provided for information only, and product specifications are subject to change without notice. Any change will be reflected in future publishings.

# Table of Contents

---

<b>1.0</b>	<b>Introduction .....</b>	<b>1</b>
1.1	Purpose.....	1
1.2	Overview.....	1
<hr/>		
<b>2.0</b>	<b>Command Interface.....</b>	<b>2</b>
2.1	Constructor.....	2
2.2	Functions.....	2
2.2.1	General Functions .....	2
2.2.2	Controller Command Functions .....	3
<hr/>		
<b>3.0</b>	<b>Python Example.....</b>	<b>15</b>
	<b>Service Form .....</b>	<b>17</b>





# Agilis-D Controller with Strain Gages Feedback CONEX-AGAP

---

## 1.0 Introduction

---

### 1.1 Purpose

The purpose of this document is to provide the method syntax of each command to communicate with the CONEX-AGAP device.

### 1.2 Overview

The Command Interface is the wrapper class that maintains a list of CONEX-AGAP instruments. It exposes methods to communicate with any CONEX-AGAP device.

These commands work in both synchronous and asynchronous mode in the NStruct environment or not. The communication is based on the NStruct server.

---

#### NOTE

Each function name is defined with the command code “AA”.

For each command function, refer to the CONEX-AGAP programmer’s manual.

---

## 2.0 Command Interface

---

### 2.1 Constructor

CONEXAGAP()

The constructor is used to create an instance of the CONEX-AGAP device.

### 2.2 Functions

#### 2.2.1 General Functions

##### 2.2.1.1 RegisterComponent

###### Syntax

int RegisterComponent(string instrumentKey)

instrumentKey: Instrument key

return: componentID

###### Description

This function allows registering the device to the server. A component ID is returned. If the registering failed, the returned component ID is zero.

---

###### NOTE

**The component ID is mandatory to use all commands from the CommandInterface.**

---

##### 2.2.1.2 UnregisterComponent

###### Syntax

void UnregisterComponent(int componentID)

###### Description

This function allows unregistering the device to the server with the component ID.

##### 2.2.1.3 LockInstrument

###### Syntax

int LockInstrument(int componentID, int stage, ref string response)

###### Description

This function allows locking the device communication to not share the communication.

##### 2.2.1.4 UnlockInstrument

###### Syntax

int UnlockInstrument(int componentID, int stage, ref string response)

**Description**

This function allows unlocking the device communication to share the communication.

**2.2.1.5 GetDevices****Description**

This function is used to get the list of the connected devices

**Syntax**

```
string[] GetDevices()
```

Return: string array to get the list of devices

**2.2.1.6 WriteToInstrument****Syntax**

```
int WriteToInstrument(int componentID, string command, ref string response, int stage)
```

componentID: Instrument ID

command: Instrument command

response: Response of the command

stage: Instrument Stage

Return: Communication error code

**Description**

This Overridden function Queries or writes the command given by the user to the instrument.

**2.2.2 Controller Command Functions****2.2.2.1 ID\_Get****Syntax**

```
int ID_Get(int componentID, int controllerAddress, out string outStageIdentifier, out string errString)
```

componentID: Instrument ID

controllerAddress: Address of Controller

outStageIdentifier: outStageIdentifier

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous ID Get command which is used to Get stage identifier.

**2.2.2.2 ID\_Set****Syntax**

int ID\_Set(int componentID, int controllerAddress, string inStageIdentifier, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inStageIdentifier: inStageIdentifier.

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous ID Set command which is used to Set stage identifier.

**2.2.2.3 JA\_Get****Syntax**

int JA\_Get(int componentID, int controllerAddress, string Axis, out double outJogVelocity, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outJogVelocity: outJogVelocity

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous JA Get command which is used to Move jogging.

**2.2.2.4 JA\_Set****Syntax**

int JA\_Set(int componentID, int controllerAddress, string Axis, double inJogVelocity, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inJogVelocity: inJogVelocity.

errString: The failure reason

Return: 0 in success and -1 on failure



**Description**

This function is used to process synchronous JA Set command which is used to Move jogging.

**2.2.2.5 MM\_Get****Syntax**

int MM\_Get(int componentID, int controllerAddress, out string outState, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outState: outState

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous MM Get command which is used to Leave DISABLE state.

**2.2.2.6 MM\_Set****Syntax**

int MM\_Set(int componentID, int controllerAddress, int inState, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inState: inState.

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous MM Set command which is used to Leave DISABLE state.

**2.2.2.7 PA\_Get****Syntax**

int PA\_Get(int componentID, int controllerAddress, string Axis, out double outTarget, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outTarget: outTarget

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PA Get command which is used to Move absolute.

### 2.2.2.8 PA\_Set

#### **Syntax**

int PA\_Set(int componentID, int controllerAddress, string Axis, double inTarget, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inTarget: inTarget.

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous PA Set command which is used to Move absolute.

### 2.2.2.9 PR\_Get

#### **Syntax**

int PR\_Get(int componentID, int controllerAddress, string Axis, out double outStep, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outStep: outStep

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous PR Get command which is used to Move relative.

### 2.2.2.10 PR\_Set

#### **Syntax**

int PR\_Set(int componentID, int controllerAddress, string Axis, double inStep, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inStep: inStep.

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous PR Set command which is used to Move relative.

**2.2.2.11 PW\_Get****Syntax**

int PW\_Get(int componentID, int controllerAddress, out int outState, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outState: outState

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PW Get command which is used to Enter/Leave CONFIGURATION state.

**2.2.2.12 PW\_Set****Syntax**

int PW\_Set(int componentID, int controllerAddress, int inState, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inState: inState.

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous PW Set command which is used to Enter/Leave CONFIGURATION state.

---

**NOTE**

**The PW command is limited to 100 writes. Unit failure due to excessive use of the PW command is not covered by warranty.**

**The PW command is used to change the configuration parameters that are stored in memory, and not parameters that are needed to be changed on the fly.**

---

**2.2.2.13 RS****Syntax**

int RS(int componentID, int controllerAddress, out string errString)

clientID: Instrument ID

controllerAddress: controllerAddress identifying the Address of Controller

errString: The failure reason

Return: 0 in success and -1 on failure

**Description**

This function is used to process synchronous RS Set command which is used to Reset controller.

#### 2.2.2.14 RS485

##### Syntax

int RS485(int componentID, int controllerAddress, out string errString)

clientID: Instrument ID

controllerAddress: controllerAddress identifying the Address of Controller

errString: The failure reason

Return: 0 in success and -1 on failure

##### Description

This function is used to process synchronous RS## Set command which is used to Reset controller's address to 1.

#### 2.2.2.15 SA\_Get

##### Syntax

int SA\_Get(int componentID, int controllerAddress, out int outAdress, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outAdress: outAdress

errString: The failure reason

Return: 0 in success and -1 on failure

##### Description

This function is used to process synchronous SA Get command which is used to Get controller's RS-485 address.

#### 2.2.2.16 SA\_Set

##### Syntax

int SA\_Set(int componentID, int controllerAddress, int inAdress, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inAdress: inAdress.

errString: The failure reason

Return: 0 in success and -1 on failure

##### Description

This function is used to process synchronous SA Set command which is used to Set controller's RS-485 address.

### 2.2.2.17 SL\_Get

#### **Syntax**

int SL\_Get(int componentID, int controllerAddress, string Axis, out double outLimit, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outLimit: outLimit

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous SL Get command which is used to Get negative software limit.

### 2.2.2.18 SL\_Set

#### **Syntax**

int SL\_Set(int componentID, int controllerAddress, string Axis, double inLimit, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inLimit: inLimit.

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous SL Set command which is used to Set negative software limit.

### 2.2.2.19 SR\_Get

#### **Syntax**

int SR\_Get(int componentID, int controllerAddress, string Axis, out double outLimit, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outLimit: outLimit

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous SR Get command which is used to Get positive software limit.

### 2.2.2.20 SR\_Set

#### **Syntax**

int SR\_Set(int componentID, int controllerAddress, string Axis, double inLimit, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inLimit: inLimit.

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous SR Set command which is used to Set positive software limit.

### 2.2.2.21 ST

#### **Syntax**

int ST(int componentID, int controllerAddress, out string errString)

clientID: Instrument ID

controllerAddress: controllerAddress identifying the Address of Controller

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous ST Set command which is used to Stop motion

### 2.2.2.22 TB

#### **Syntax**

int TB(int componentID, int controllerAddress, string inErrorCode, out string outErrorCode, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inErrorCode: inErrorCode.

outErrorCode: outErrorCode

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TB Get command which is used to Get command error string.

### 2.2.2.23 TE

#### **Syntax**

int TE(int componentID, int controllerAddress, out string outLastCommandError, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outLastCommandError: outLastCommandError

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TE Get command which is used to Get last command error.

### 2.2.2.24 TH

#### **Syntax**

int TH(int componentID, int controllerAddress, string Axis, out double outPosition, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outPosition: outPosition

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TH Get command which is used to Get target position.

### 2.2.2.25 TP

#### **Syntax**

int TP(int componentID, int controllerAddress, string Axis, out double outPosition, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outPosition: outPosition

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TP Get command which is used to Get current position.

### 2.2.2.26 TS

#### **Syntax**

int TS(int componentID, int controllerAddress, out string errorCode, out string controllerState, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

errorCode: errorCode

controllerState: controllerState

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous TS Get command which is used to Get positioner error and controller state.

### 2.2.2.27 VE

#### **Syntax**

int VE(int componentID, int controllerAddress, out string outInformation, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outInformation: outInformation

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous VE Get command which is used to Get controller revision information.

### 2.2.2.28 XU\_Get

#### **Syntax**

int XU\_Get(int componentID, int controllerAddress, string Axis, out double outStepSize, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outStepSize: outStepSize

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous XU Get command which is used to Get step size for STEPPING OL state.



### 2.2.2.29 XU\_Set

#### Syntax

int XU\_Set(int componentID, int controllerAddress, string Axis, double inStepSize, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inStepSize: inStepSize.

errString: The failure reason

Return: 0 in success and -1 on failure

#### Description

This function is used to process synchronous XU Set command which is used to Set step size for STEPPING OL state.

### 2.2.2.30 XR\_Get

#### Syntax

int XR\_Get(int componentID, int controllerAddress, string Axis, out double outStep, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

outStep: outStep

errString: The failure reason

Return: 0 in success and -1 on failure

#### Description

This function is used to process synchronous XR Get command which is used to Move stepping.

### 2.2.2.31 XR\_Set

#### Syntax

int XR\_Set(int componentID, int controllerAddress, string Axis, double inStep, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

inStep: inStep.

errString: The failure reason

Return: 0 in success and -1 on failure

#### Description

This function is used to process synchronous XR Set command which is used to Move stepping.

### 2.2.2.32 ZT

#### **Syntax**

int ZT(int componentID, int controllerAddress, out List<string> Parameters, out string errString)

componentID: Instrument ID

controllerAddress: Address of Controller

Parameters: Parameters

errString: The failure reason

Return: 0 in success and -1 on failure

#### **Description**

This function is used to process synchronous ZT Get command which is used to Get all controller parameters.

### 3.0 Python Example

```

=====
#Initialization Start
#The script within Initialization Start and Initialization End is needed for properly
# initializing IOPortClientLib and Command Interface for CONEX-AGAP
instrument.
#The user should copy this code as is and specify correct paths here.
import sys
import time

#IOPortClientLib and Command Interface DLL can be found here.
print "Adding location of IOPortClientLib.dll &
Newport.CONEXAGAP.CommandInterface.dll to sys.path"
sys.path.append(r'C:\Program Files\Newport\Instrument
Manager\NStruct\Instruments\CONEX-AGAP\Bin')

# The CLR module provide functions for interacting with the underlying
# .NET runtime
import clr

# Add reference to assembly and import names from namespace
clr.AddReferenceToFile("Newport.CONEXAGAP.CommandInterface.dll")
from CommandInterface import *

import System
=====

# Instrument Initialization
# The key should have double slashes since
# (one of them is escape character)
instrument="CONEX-AGAP (A6TOO38G)"
print 'Instrument Key=>', instrument

# create AGAP instance
AGAP = CONEXAGAP()

# register to server, componentID needs to be used in all commands
componentID = AGAP.RegisterComponent(instrument);
print 'componentID=>', componentID

# variable setting
address = 1
axis = "U"

# Get controller revision information
result, response, errString = AGAP.VE(componentID,address)
if result == 0 :
    print 'controller revision=>', response
else:
    print 'Error=>',errString

# Get positive software limit
result, response, errString = AGAP.SR_Get(componentID,address,axis)
if result == 0 :
    print 'positive software limit=>', response
else:
    print 'Error=>',errString

```

```
# Get negative software limit
result, response, errString = AGAP.SL_Get(componentID,address,axis)
if result == 0 :
    print 'negative software limit=>', response
else:
    print 'Error=>',errString

# Get current position
result, response, errString = AGAP.TP(componentID,address,axis)
if result == 0 :
    print 'position=>', response
else:
    print 'Error=>',errString

print 'End of script'

# unregister server
AGAP.UnregisterComponent(componentID);
```





Visit Newport Online at:  
[www.newport.com](http://www.newport.com)

**North America & Asia**

Newport Corporation  
1791 Deere Ave.  
Irvine, CA 92606, USA

**Sales**

Tel.: (800) 222-6440  
e-mail: [sales@newport.com](mailto:sales@newport.com)

**Technical Support**

Tel.: (800) 222-6440  
e-mail: [tech@newport.com](mailto:tech@newport.com)

**Service, RMAs & Returns**

Tel.: (800) 222-6440  
e-mail: [service@newport.com](mailto:service@newport.com)

**Europe**

MICRO-CONTROLE Spectra-Physics S.A.S  
9, rue du Bois Sauvage  
91055 Évry CEDEX  
France

**Sales**

Tel.: +33 (0)1.60.91.68.68  
e-mail: [france@newport.com](mailto:france@newport.com)

**Technical Support**

e-mail: [tech\\_europe@newport.com](mailto:tech_europe@newport.com)

**Service & Returns**

Tel.: +33 (0)2.38.40.51.55

