


XPS

Universal High-Performance Motion Controller/Driver



Newport®

Experience | Solutions

Stages Database History

Precision Motion – **Guaranteed™**

Preface

Confidentiality & Proprietary Rights

Reservation of Title

The Newport Programs and all materials furnished or produced in connection with them ("Related Materials") contain trade secrets of Newport and are for use only in the manner expressly permitted. Newport claims and reserves all rights and benefits afforded under law in the Programs provided by Newport Corporation.

Newport shall retain full ownership of Intellectual Property Rights in and to all development, process, align or assembly technologies developed and other derivative work that may be developed by Newport. Customer shall not challenge, or cause any third party to challenge, the rights of Newport.

Preservation of Secrecy and Confidentiality and Restrictions to Access

Customer shall protect the Newport Programs and Related Materials as trade secrets of Newport, and shall devote its best efforts to ensure that all its personnel protect the Newport Programs as trade secrets of Newport Corporation. Customer shall not at any time disclose Newport's trade secrets to any other person, firm, organization, or employee that does not need (consistent with Customer's right of use hereunder) to obtain access to the Newport Programs and Related Materials. These restrictions shall not apply to information (1) generally known to the public or obtainable from public sources; (2) readily apparent from the keyboard operations, visual display, or output reports of the Programs; (3) previously in the possession of Customer or subsequently developed or acquired without reliance on the Newport Programs; or (4) approved by Newport for release without restriction.

©2016 Newport Corporation
1791 Deere Ave.
Irvine, CA 92606, USA
(949) 863-3144

XPS

Universal High-Performance Motion Controller/Driver

Versions 3.0.x

XPS Stages Database Version: V3.0.6

Version history: September, 19th 2016 (ECO 4552).

- Add new stages on IDL family :
 - Standard configurations with XPS-EDBL for:
IDL225-400LM, IDL225-1200LM, IDL280-400LM, IDL280-1200LM, IDL560-450LM
 - LMI configurations with XPS-EDBL for:
IDL225-400LM, IDL225-1200LM, IDL280-400LM, IDL280-1200LM, IDL560-450LM
- Modification on RGV160BL configurations for XPS-EDBL :
 - No SmartStageName

XPS Stages Database Version: V3.0.5

Version history: December, 11th 2015 (ECO 4490).

- Add new parameters so the Database is now compatible with the XPS-RL controller.
- Remove obsolete sections (PreFeedForward filter et PostFeedForward filter).
- Add a checksum to verify the integrity of the 'XPS-StageDataBase' file.
- Add a new configuration : DUMMY@DUMMY_STAGE@NO_DRIVER
- Modification of ONE-XY family :
 - HomeSearch.MechanicalZeroAndIndexHomeSearch is replaced by HomeSearch.MinusEndOfRunAndIndexHomeSearch
 - HomeSearchMaximumAcceleration was 100 and is set to 250
 - DriverCutOffFrequency was 300 and is set to 200
 - **ONE-XY@ONE-XY300-X :**
 - ScalingAcceleration was 1605 and is set to 1497
 - AccelerationLimit was 1494 and is set to 1393

- ONE-XY@ONE-XY300-Y :
- ScalingAcceleration was 3351 and is set to 3125
- AccelerationLimit was 3118 and is set to 2908

Versions 2.8.x

XPS Stages Database Version: **V2.8.12**

Version history: December, 11th 2015 (ECO4429).

- Modification of MTN@MTN200PP@XPS-DRV01: SmartStageName was MTN100PP and is set to MTN200PP (Wrong SmartStageName)

XPS Stages Database Version: **V2.8.11**

Version history: November, 2nd 2015 (ECO4379).

- Add ONE-XY family: ONE-XY60, ONE-XY100, ONE-XY200, ONE-XY300 configurations with DRV02and DRV02-LMI.

XPS Stages Database Version: **V2.8.10**

Version history: October, 15th 2015 (ECO4370).

- Modification on configurations XPS-DRV03 for MTN100CC, MTN200CC, MTN300CC:
 - ScalingCurrent was 3 and is set to 5
 - CurrentLimit was 3 and is set to 4.9.
- Add IDL family:
 - Standard configurations with XPS-EDBL for:
 - IDL165-200LM, IDL225-300LM IDL225-600LM, IDL280-300LM, IDL280-600LM, IDL560-600LM
 - LMI configurations with XPS-EDBL for:
 - IDL165-200LM, IDL225-300LM IDL225-600LM, IDL280-300LM, IDL280-600LM, IDL560-600LM
 - Standard configuration with XPS-DRV02 for IDL165-200LM@XPS-DRV02
 - LMI configuration with XPS-DRV02 for IDL165-200LM@XPS-DRV02

XPS Stages Database Version: **V2.8.9**

Version history: June, 25th 2015 (ECO4329).

- Modification on configurations XPS-DRV01 for MTN100CC, MTN200CC, MTN300CC:
CurrentLimit was 3.0A and is set to 2.9A.
- Modification on configurations XPS-DRV01 for MTN100PPV6, MTN200PPV6 and MTN300PPV6:
PeakCurrentPerPhase was 1.25A and is set to 1.8A.

XPS Stages Database Version: V2.8.8**Version history: April, 23th 2015 (ECO4306).**

- Add MTN100PPV6, MTN200PPV6 and MTN300PPV6 configurations for XPS-DRV01

XPS Stages Database Version: V2.8.7**Version history: February, 24th 2015 (ECO4268).**

Modifications

- Modification on configurations XPS-DRV01 and XPS-DRV03 for MTN100CC, MTN200CC, MTN300CC:
VoltageLimit was 38.4V and is set to 45V.

XPS Stages Database Version: V2.8.6**Version history: January, 12th 2015 (ECO3478 and ECO4266).**

Modifications

- Add configurations: RV@RV120HAHL@XPS-DRV03, RV@RV120HAHL-F@XPS-DRV03

XPS Stages Database Version: V2.8.5**Version history: January 2015 (ECO4237).**

Modifications

- Modification on RVS80CC, UTS50CC, UTS100CC, UTS150CC, FMS100CC, FMS200CC, FMS300CC configurations:
DerivativeFilterCutOffFrequency was 4000 and is set to 1000 Hz and Kd becomes 0.8.
- Add configurations on 'Piezo' Family: NPA25-D, NPA50-D, NPA100-D, NPM140-D, NP0140-D, NPO250-D, NPXY100-D-X, NPXY100-D-Y, NPXY100-D-Z, NPXY200-D-X, NPXY200-D-Y, NPXYZ100-D-X, NPXYZ100-D-Y, NPXYZ100-D-Z

XPS Stages Database Version: V2.8.4**Version history: August 2014 (ECO4187).**

Modifications

- Modification on NPX200-D / NPX400-D / PSM2-D-X / PSM2-D-Y / PSM2-D-Z configurations:
DriverLowpassFilter was 0 and is set to 2000
- Replace configurations remove on the V2.8.2

XPS Stages Database Version: V2.8.3**Version history: August 2014 (ECO4174).**

Modifications

- Add IMSxx-LM and IMSxx-LM-SA configurations for XPS-EDBL
- Add RGV100HL configurations for XPS-EDBL
- Add RGV160BL configurations for XPS-EDBL

XPS Stages Database Version: V2.8.2

Version history: July 2014 (ECO4171).

Modifications

- Modification on RGV100HL configurations: DerivativeFilterCutOffFrequency was 4000 and is set to 500
- Add MTN100CC, MTN200CC and MTN300CC configurations for XPS-DRV01
- Add MTN100CC, MTN200CC and MTN300CC configurations for XPS-DRV03
- Add MTN100PP, MTN200PP and MTN300PP configurations for XPS-DRV01

XPS Stages Database Version: V2.8.1

Version history: Septembre 2013 (ECO4061).

Modifications

- Delete TRB6PPD and TRB12PPD configurations
- Add TRB6PP, TRB12PP, TRB25PP configurations

XPS Stages Database Version: V2.8.0

Version history: Septembre 2013 (ECO4056).

Modifications

- Update section's names.
- Remove obsolete configurations.
- Rename "StandardEOREncoderPlug" by "StandardLimitAndHomeEncoderPlug".

Versions 2.7.x

XPS Stages Database Version: V2.7.7

Version history: July 2013 (ECO4043).

Modifications

- Modification on URB100CC configurations (DRV01, DRV01-Spindle, DRV03 and DRV03-Spindle):
Maximum allowed motor voltage 14.56V to 20V
- Modification on URS50BCC and URS50BPP for spindle configuration:
'ServitudesType' become 'Spindle'

XPS Stages Database Version: V2.7.6

Version history: May 2013 (ECO 4015).

All stages

- Add NPM140-SG driven by DRV1.
- Add TRB family: TRB6CC, TRB6PPD, TRB12CC, TRB12PPD, TRB25CC, TRB25PPD

Modifications

- Modification on 'EncoderZMPlug = Driver' for RGV100HL configurations.
- Modification on 'EncoderIndexOffset =0' for FMSxxxPPHA configurations.
- Modification of 'Initialization acceleration level' = 25% on ILSxxLM configuration

XPS Stages Database Version: [V2.7.5](#)**Version history: April 2013 (ECO 3999).**

Not released.

XPS Stages Database Version: [V2.7.4](#)**Version history: July 2012 (ECO 3848).**

All stages

- Add RGV100HL configurations driven by DRV02.
- Add FMS100PPHA, FMS200PPHA, FMS300PPHA configurations driven by DRV01.

Modifications

- Modification of « DerivativeFilterCutOffFrequency » value for all configurations superior at 4000 and became equal to 4000.

XPS Stages Database Version: [V2.7.3](#)**Version history: May, 2012 (ECO 3811).**

All stages

- Add NPA50SG-D, NPXY100SG-D-X, NPXY100SG-D-Y, NPXY200SG-D-X, NPXY200SG-D-Y configurations with DRV1.

XPS Stages Database Version: [V2.7.2](#)**Version history: March, 2012 (ECO 3783).**

All stages

- Add NPX200-D with DRV1.

Modifications

- Modification of PSM2-D-X, PSM2-D-Y, and PSM2-D-Z configurations with DRV1.

XPS Stages Database Version: [V2.7.1](#)**Version history: January, 2012 (ECO 3756).**

All stages

- Add NPA25SG-D AND NPX400SG-D with DRVP1.
- Add URB100CC configuration for DRV01 and DRV01-Spindle

Modifications

- Modification of the IntegrationTime for all RVxxHAHLT, RVxxHAHLT-F, RVxxHAT and RVxxHAT-F configurations driven by DRV01, DRV03 and DRVM3: new value equal to 3s
- Modification of the DeadBandThreshold for all RVxxHAHLT, RVxxHAHLT-F, RVxxHAT and RVxxHAT-F configurations driven by DRV01, DRV03 and DRVM3: new value equal to 0.

XPS Stages Database Version: **V2.7.0**

Version history: December, 2011 (ECO 3739).

All stages

- Add PSM2-D-XYZ with DRVP1

Modifications for QNX firm

- **EncoderIndexOffset:** This new option is used for Interpolated and AquadB Encoders
- **EncoderHardInterpolatorErrorCheck:** This new option is used for Interpolated Encoder
- **HomingSensorOffset:** This new option is used for every types of Home search processes, unless CurrentPositionAsHome

Versions 2.6.x

XPS Stages Database Version: **V2.6.9**

Version history: November 8, 2011 (ECO 3716).

All stages

- Add ZVR-PP-T, ZVR-PP-Z, ZVR-PC-T, ZVR-PC-Z

XPS Stages Database Version: **V2.6.8**

Version history: October 20, 2011 (ECO 3704).

All stages

- Add IMS800LM-SA, IMS1000LM-SA, IMS1200LM-SA, and IMS800LM, IMS1000LM, IMS1200LM with DRV02 and DRV02-LMI
- Add NPX200SG-D with DRVP1

Modifications

- Modifications on TRA25CC with DRV01: Maximum allowed motor current 1.5A to 0.15 A

XPS Stages Database Version: **V2.6.7**

Version history: June 17, 2011 (ECO 3618).

All stages

- Add URS50BPPV6 with DRV01
- Add LTAHLPPV6 & LTAHSPPV6

Modifications

- Modifications on FMS100PP, FMS200PP, FMS300PP: No encoder and No servo loop with position output
- Modifications on URB100CC: End of Run = +/- 165°

XPS Stages Database Version: V2.6.6**Version history: March 29, 2011.**

All stages

- Add FMS100CC, FMS200CC, FMS300CC and FMS100PP, FMS200PP, FMS300PP with DRV01

Modifications

- Modifications on GTS70 et GTS150:
 - Maximum allowed motor voltage 25V to 36V
 - Maximum allowed motor current 2A to 1.6A

XPS Stages Database Version: V2.6.5**Version history: January 24, 2011.**

Modifications

- Modify the following error of NPXYZ100SG-D-X, NPXYZ100SG-D-Y, NPXYZ100SG-D-Z
- Add NPA100SG-D, NPO250SG-D, NPX400-D

XPS Stages Database Version: V2.6.4**Version history: December 17, 2010.**

Modifications

- Add NPXYZ100SG-D-X, NPXYZ100SG-D-Y, NPXYZ100SG-D-Z

XPS Stages Database Version: V2.6.3**Version history: November 25, 2010.**

Modifications

- ILS50HA, ILS100HA, ILS150HA, ILS200HA, ILS250HA with DRV03 in closed loop
- RV350CC-F with DRV03 in closed loop
- RV350CC with DRV03 in closed loop
- RV350HAHLT with DRVM3 in closed loop

XPS Stages Database Version: V2.6.2**Version history: November 23, 2010.**

All stages

- Add ILS50HA, ILS100HA, ILS150HA, ILS200HA, ILS250HA with DRV01 and DRV03

Modifications

- Change ILS-LM with ILSLM

Deleted configurations

XPS Stages Database Version: V2.6.1**Version history: September 1, 2010.**

All stages

- Add ILS100-LM, ILS200-LM, ILS300-LM
- Add RV120PE-FBV6 RV120PEBV6 (and spindle mode)

Modifications

- Modify the velocity servo loop cut off frequency of RV350HAHLT/RV350HAHLT-F driven by DRV03
old value 100 Hz to new value 10 Hz
- Modify PID parameters of RV350CC/RV350CC-F driven by DRV03
KP = 100, KI = 2000, KD=1, Fc=1000

Deleted configurations

XPS Stages Database Version: V2.6.0**Version history: August 19, 2010.**

Change of name because of V2.6 firmware release.

No change of the contents.

Versions 2.5.x**XPS Stages Database Version: V2.5.4****Version history: June 25, 2010.**

All stages

- Add URS50BCC and URS50BPP (DRV01 and DRV01@Spindle configurations)

Modifications

Deleted configurations

XPS Stages Database Version: V2.5.3**Version history: February 25, 2010.**

All stages

- Add DRV03 configuration for TRA12CC
- Add configuration DRV01 for NSA12
- Add configuration DRV03 for LTA-HXP100
- Add configurations DRV03 for URB100CC (Spindle and with EoR)

Modifications

- ILS-CCHA: Modifications of PID parameters ($K_P = 100$, $K_I = 5000$, $K_D = 0.1$) and $DerivativeFilterCutOffFrequency = 500$
- VP-25AA, VP-25XA, VP-25XL, VP-5ZA: Modifications of $ScalingVelocity$ and $VelocityLimit$ equal to $V_{max} + 50\%$ for DRV03 configurations

Deleted configurations

XPS Stages Database Version: [V2.5.1 & V2.5.2](#)

Not applicable

XPS Stages Database Version: [V2.5.0](#)

Version history: June 29, 2009.

Change of name because of V2.5 firmware release.

No change of the contents.

Versions 2.4.x

XPS Stages Database Version: [V2.4.3](#)

Version history: ECO released 04, 2009

All stages

- Add DRV03 configurations missing in V2.4.2

Modifications

- TRA6CC and TRA25CC DRV03: $ClosedLoopStatus = Closed$

Deleted configurations

- RV120HAT Spindle
- TRA12CC DRV01

XPS Stages Database Version: V2.4.2**Version history:**

All stages

- PR50CC: HomeSearchTimeOut from 400003 s to 40s
- In “PID with acceleration output” configuration: new parameter “KFeedForwardJerk=0”
- DRV03 Configuration for RV120HAHLT and RV120HAHLT-F
- Modifications of PID parameters ($K_p=500$, $K_i=50$) and Servo loop dead band threshold = $2 * \text{resolution}$ for:
 - RV120/160/240/350/HAT
 - RV120/160/240/350/HAT-F
 - RV160/240/350/HAHLT
 - RV160/240/350/HAHLT-F

New stages

- TRA6PPV6, TRA12PPV6 and TRA25PPV6
- Stages with new UE41PPV6 motor (Half becomes full, and Current / 2)
 - MTM100/150/200/250/PE1BV6
 - MTM100/150/200/250/PE.1BV6
 - BGM120/160/PEBV6
 - UZM160PE.05BV6

XPS Stages Database Version: V2.4.1**Version history: Web released 11, 2008**

All stages

- TRA6CC and TRA25CC following error changed from 1 to 0.01 mm
- PR50CC and SR50CC encoder resolution changed in spindle mode to be a submultiple of 360°
- PR50PP and SR50PP displacement per full step changed in spindle mode to be a submultiple of 360°

New stages

- None

XPS Stages Database Version: V2.4.0**Version history: ECO released 07, 2008**

All stages

- MFA-CC and MFA-CCV6
 - travel range changed to 0 – 25 mm
 - PID parameters optimized

New stages

- MFA-PPD with XPS-DRV01 driver board added
- TRA6CC and TR25CC with XPS-DRV01 and XPS-DRV03 driver boards added

- TRA6PPD and TRA25PPD with XPS-DRV01 driver board added
- TRA6PPV6 and TRA25PPV6 with XPS-DRV01 driver board added

Versions 2.3.x

XPS Stages Database Version: V2.3.0

Version history: ECO released 05, 2008

All stages

- Missing closed loop status parameter for IMSxxV stages corrected

New stages

- None

Versions 2.2.x

XPS Stages Database Version: V2.2.0

Version history: ECO released March, 2008

All stages

- Minus limit changed for CMAxx from 0 to -0.0001 according to the new constraint of the V2.2.0 firmware
- Minus limit changed for LTAXx from 0 to -0.0001 according to the new constraint of the V2.2.0 firmware
- Minus limit changed for UTMxx/-1 from 0 to -0.0001 according to the new constraint of the V2.2.0 firmware

New stages

- None

Versions 2.1.x

XPS Stages Database Version: V2.1.3

Version history: Web released March, 2008

All stages

- Missing closed loop status parameter for RV stages with XPS-DRV03 driver board corrected
- KP and GKP parameters changed for IMSxx-V from 1000 to 800 and from 6 to 7.5
- XMSxx and XMLxx home search timeout increased

- Current servo loop bandwidth of RVxxT stages with XPS-DRV03 driver board changed from 750 to 500 Hz
- Velocity servo loop bandwidth of RVxxT stages with XPS-DRV03 driver board changed from 150 to 100 Hz

New stages

- None

XPS Stages Database Version: V2.1.2

Version history: ECO released November, 2007

All stages

- Default value of KForm parameter changed to 0
- XMS and XML home search velocity changed to 100 mm/s and time out updated
- BGS80PP limits corrected
- UTSxxCC minimum jerk time and maximum jerk time changed from 0.001 and 0.01 to 0.005 and 0.05
- LTA and BGS backlash changed to 0 unit
- ScalingCurrent corrected to 5 A for XPS-DRV03 driver board configuration
- CurrentLimit corrected for XPS-DRVMx driver board configuration

New stages

- UZS80CC and UZS80PP with XPS-DRV01 and XPS-DRV03 (CC only) driver boards added
- URSxxBCC, URSxxBPP and URS75BPPV6 with XPS-DRV01 and XPS-DRV03 (CC only) driver boards added

XPS Stages Database Version: V2.1.1

Version history: ECO released July, 2007

All stages

- VP-5ZA travel changed to +/- 2.4 mm
- RV120CCHL, RV160CCHL, RV240CCHL, RV350CCHL servo loop gains corrected
- RV120CC, RV160CC, RV240CC, RV350CC with XPS-DRV03 driver board servo loop gains corrected
- BGS50CC encoder resolution rounding corrected
- BGS50PP displacement per full step rounding corrected
- UTSCC maximum allowed voltage corrected
- RGV100BL with XPS-DRV02 driver board configuration name corrected
- RGV100 current cut-off frequency changed to 500 Hz

New stages

- RV350HAT and RV350HAT-F with XPS-DRV03 driver board added

XPS Stages Database Version: V2.1.0**Version history: ECO released February, 2007**

All stages

- LTAPPV6 home search timeout increased to 400 s
- PR50CC encoder resolution corrected
- XMS, XML and IMS-LM with AnalogSin120AccelerationLMI interface added
- InitializationAcceleration parameters replaced by InitializationAccelerationLevel in %: need a manual update by the user for current stages.ini configurations
- MappingFileName parameter removed for NoEncoder configuration

New stages

- UTS50CC, UTS100CC, UTS150CC, UTS50PP, UTS100PP, UTS150PP, UTS50PPV6, UTS100PPV6 et UTS150PPV6 with XPS-DRV01 and XPS-DRV03 (CC only) driver boards added
- RVS80CC et RVS80PP with XPS-DRV01 and XPS-DRV03 (CC only) driver boards added
- GTS30V with XPS-DRV01 and XPS-DRV03 driver boards added
- RGV100BL with XPS-DRV02 driver board added
- BGS50CC, BGS50PP, BGS80CC, BGS80PP with XPS-DRV01 and XPS-DRV03 (CC only) driver boards added
- URS75PPV6 with XPS-DRV01 driver board added

Versions 2.0.x**XPS Stages Database Version: V2.0.0****Version history: Released November, 2006**

All stages

- GTS missing parameters added
- M- configurations removed to reduce stage data base size
- XPS-DRV03 configuration added for all ILSCCHA stages

Versions 1.6.x**XPS Stages Database Version: V1.6.3****Version history: Released October, 2006**

All stages

- DriverCurrentCutOffFrequency parameter name change to DriverCutOffFrequency for XPS-DRV02 driver board
- PositionerMappingLineNumber and PositionerMaxPositionError parameters added in case of mapping use but empty

XPS Stages Database Version: V1.6.2**Version history: Released October, 2006**

All stages

- (M-)UTM stepper /1 and /-1 configurations corrected
- CMA-25CCCL and CMA-12CCCL pulse width frequency changed to 200 kHz
- GTS70 and GTS150 parameters corrected
- XPS-DRVM5 and XPS-DRV03 configurations added to (M-)VP-25AA, (M-)VP-25XA, (M-)VP-5ZA
- Travel limits corrected for LTA actuators
- first release with the stage database tool manager

New stages

- LTAPPV6 added

XPS Stages Database Version: V1.6.1**Version history: Released July, 2006**

All stages

- Modify Kd for CMA12CCCL

XPS Stages Database Version: V1.6.0**Version history: Released June, 2006**

All stages

- Add XPS-DRV03 configurations in parallel to XPS-DRVMx configurations
- Add "spindle" configurations for all stages in rotation not for BGM and RV-HA
- modification PID paramètres of CMA12CCCL and CMA25CCCL
- correction PID parameter units of RGV

New stages

- GTS75 and GTS150 stages

Versions 1.5.x**XPS Stages Database Version: V1.5.0****Version history: Released January, 2006**

All stages

- Modifications of LTA
- Modifications of RV240 and RV350

Versions 1.4.x

XPS Stages Database Version: [V1.4.3](#)

Version history: Released November, 2005

All stages

- Correction of a bug concerning the travel of all DRVM driven stages

XPS Stages Database Version: [V1.4.2](#)

Version history: Released November, 2005

All stages

- Modification of DriverErrorAmplifierGain from 9 to 5 on all MTM???CC?? stages

New stages

- RV stages with -F extension
- M-IMS100V and M-IMS300V

XPS Stages Database Version: [V1.4.1](#)

Version history: Version July, 2005

All stages

- Usage of variable PID
- Short integration time: 1sec in place of 1e+99sec
- jerktime = 0.5sec in place of 0.05
- HomeSearchMaximumAcceleration = MaximumAcceleration/4 in place of MaximumAcceleration/2

New stages

- IMS100V
- IMS300V

XPS Stages Database Version: [V1.4.0](#)


Version history: Released May, 2005

All stages

- None

New stages

- M-MTM100PE.1 encoder resolution corrected
- IMS LM added
- XML and XMS added



Newport®
Experience | Solutions

Visit Newport Online at:
www.newport.com

North America & Asia

Newport Corporation
1791 Deere Ave.
Irvine, CA 92606, USA

Sales

Tel.: (800) 222-6440
e-mail: sales@newport.com

Technical Support

Tel.: (800) 222-6440
e-mail: tech@newport.com

Service, RMAs & Returns

Tel.: (800) 222-6440
e-mail: service@newport.com

Europe

MICRO-CONTROLE Spectra-Physics S.A.S
9, rue du Bois Sauvage
91055 Évry CEDEX
France

Sales

Tel.: +33 (0)1.60.91.68.68
e-mail: france@newport.com

Technical Support

e-mail: tech_europe@newport.com

Service & Returns

Tel.: +33 (0)2.38.40.51.55