




News  
that will  
change  
the shape  
of things  
to come.



Imagine what you could do with a programmable tool that reshapes the spectral content of input light on demand, and in real time. Dynamic gain equalization in telecom amplifiers...fiber optic sensing...optical coherence tomography — the possibilities are endless.

Newport makes this all a reality with a completely new class of instruments called Programmable Spectral Processors. The first of these innovative products is the  $\lambda$ Commander™. This fiber-coupled instrument can apply any arbitrary spectral transformation to light — from simply passing selected narrow wavelengths and blocking others, to applying a complex, time varying function to the spectral shape with bandwidths up to 100 nm.



The  $\lambda$ Commander is just the start. Stay tuned for more products based on this

revolutionary technology, which will totally reshape the way you work with light. To learn more, visit

[www.newport.com/lambdacommander](http://www.newport.com/lambdacommander)  
or call 800-222-6440, ext 31354.

