



ESP302

Friendly Motion Controller/Driver



Stages Database
History

©2020 by Newport Corporation, Irvine, CA. All rights reserved.

Original instructions.

No part of this document may be reproduced or copied without the prior written approval of Newport Corporation. This document is provided for information only, and product specifications are subject to change without notice. Any change will be reflected in future publishings.



Friendly Motion Controller/Driver

Versions 4.3.x

ESP302 Stages Database Version: V4.3.15

Version history: January, 17th 2022 (ECO 5639)

- Reintroducing the configuration ZVR@ZVR-PP-Z@XPS-DRV11L

ESP302 Stages Database Version: V4.3.11

Version history: January, 17th 2022 (ECO 5510)

- Limitation of « DriverBemfCompensationGain » parameter to 1.5 for Stepper configurations below with XPS-DRV11L :
 - BGM@BGMxxxBPE@ XPS-DRV11L
 - BGM@BGMxxxBPP@XPS-DRV11L
 - BGM@BGMxxxPEBV6@XPS-DRV11L
 - BGS@BGS80PP@XPS-DRV11L
 - FMS@FMSxxxPP@XPS-DRV11L
 - FMS@FMSxxxPPHA@XPS-DRV11L
 - ILS@ILSBPP@xxxXPS-DRV11L
 - IMS@IMSxxxBPP@XPS-DRV11L
 - MTN@MTNxxxPP@XPS-DRV11L
 - MTN@MTNxxxPPV6@XPS-DRV11L
 - RV@RVxxxBPE-F@XPS-DRV11L
 - RV@RVxxxBPE@XPS-DRV11L
 - RV@RVxxxBPEV6-F@XPS-DRV11L
 - RV@RVxxxBPEV6@XPS-DRV11L
 - RV@RVxxxBPP-F@XPS-DRV11L
 - RV@RVxxxBPP@XPS-DRV11L
 - RV@RVxxxPE-FBV6@XPS-DRV11L
 - RV@RVxxxPEBV6@XPS-DRV11L
 - RVS@RVS80PP@XPS-DRV11L
 - URS@URSxxxBPP@XPS-DRV11L
 - URS@URS50CPP@XPS-DRV11L
 - URS@URS75BPPV6@XPS-DRV11L
 - UTS@UTSxxxPP@XPS-DRV11L
 - UTS@UTSxxxPPV6@XPS-DRV11L
 - UZM@UZM160BPE.05@XPS-DRV11L
 - UZM@UZM160BPP.05@XPS-DRV11L
 - UZM@UZM160PE.05BV6@XPS-DRV11L
 - UZS@UZS80PP@XPS-DRV11L

- ZVR@ZVR-PP-T@XPS-DRV11L

ESP302 Stages Database Version: **V4.3.5**

Version history: Released April 28th, 2020

- Introduction of all configurations compatible with the ESP302 controller with the driver @ XPS-DRV11L.

Your Local Representative

Fax: _____

Return authorization #: _____

Date: _____

Phone Number: _____

Fax Number: _____

Serial #: _____

Reasons of return of goods (please list any specific problems): _____

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



Visit Newport Online at:
www.newport.com

North America & Asia

Newport Corporation
1791 Deere Ave.
Irvine, CA 92606, USA

Sales

Tel.: (800) 222-6440
e-mail: sales@newport.com

Technical Support

Tel.: (800) 222-6440
e-mail: tech@newport.com

Service, RMAs & Returns

Tel.: (800) 222-6440
e-mail: service@newport.com

Europe

MICRO-CONTROLE Spectra-Physics S.A.S
9, rue du Bois Sauvage
91055 Évry CEDEX
France

Sales

Tel.: +33 (0)1.60.91.68.68
e-mail: france@newport.com

Technical Support

e-mail: tech_europe@newport.com

Service & Returns

Tel.: +33 (0)2.38.40.51.55



 Newport®

 Ophir®

 Spectra-Physics®