


# 1 mW 543.5 nm (GREEN) HELIUM NEON LASER MODEL: N-LGR-173

OUTPUT SPECIFICATIONS			
CW Power Output (mW)		> 1.0	
Wavelength (nm)		543.5	
Transverse Mode		> 90% TEM <sub>00</sub>	
Polarization		Random	
Beam Diameter at 1/e <sup>2</sup> Points		0.79 ± 5%	
Beam Divergence (mrad)		0.88 ± 5%	
Longitudinal Mode Spacing (MHz)		373	
Mode Sweeping		< 5%	
Long Term Power Drift (8 hrs)		< 2.5%	
Amplitude Noise, 30 Hz to 30 MHz (peak-to-peak)		< 2.8%	
Warmup to > 95% of Maximum Power (minutes)		< 15	
Beam Position Drift From Cold Start (mrad)		< 0.05	
Beam Position Drift After 15 Minute Warmup (mrad)		< 0.03	
ELECTRICAL SPECIFICATIONS			
Start Voltage (kVdc)		< 10	
Recommended Operating Current (mA)		6.5 ± 0.2	
Operating Voltage (VDC)		2540 ± 100	
ENVIRONMENTAL SPECIFICATIONS		OPERATING	NON-OPERATING
Temperature (°C)		-20 to +40	-40 to +80
Altitude (meters)		0 to 3000	0 to ∞
Relative Humidity (% , non-condensing)		0 to 99%	0 to 99%
Mechanical Shock (g)		< 1 for < 11 msec	< 25 for < 11 msec < 100 for < 1 msec

Specifications are subject to change without notice.



# LABORATORY HELIUM NEON LASER POWER SUPPLY

## INPUT SPECIFICATIONS

Voltage (VAC)*	115 / 230
Line Frequency (Hz)	50 to 400
Current (A)	0.2 / 0.1

## OUTPUT SPECIFICATIONS


Sustaining Voltage (VDC)	2400 to 3200
Start Voltage (kVDC)	>10
Current Setting (mA)	6.5 ± 0.2
Power (W)	< 19.1
Current Ripple (% Peak to Peak)	< 2.0
Current Ripple (% RMS)	< 0.71
Time Delay (Seconds)	3 to 7

## ENVIRONMENTAL SPECIFICATIONS

	OPERATING	NON-OPERATING
Temperature (°C)	-20 to +40	-40 to +80
Altitude (meters)	0 to 3000	0 to ∞
Relative Humidity (% , non-condensing)	0 to 99%	0 to 99%
Mechanical Shock (g)	< 50 for < 11 msec	< 50 for < 11 msec < 100 for < 1 msec

Specifications are subject to change without notice.

\*Please specify AC power cord plug type: NEMA 5-15P for 100 to 120 VAC, Europlug (CEE7/4) for 230 VAC, or British Standard (BS 1363) for 230 to 240 VAC.



Dimensions in Inches (mm)  
Reference Dimensions Only



© Newport. All Rights Reserved