


1/4-20 x 1" SPACING, 1/2" BORDER FOR VISXXXX
 M6 X 25mm SPACING, 12.5mm BORDER FOR M-VISXXXX



VIS SYSTEM	DIM A	DIM B	DIM C in.	DIM D in.
VIS3636	36.0 in	36.0 in	38.8	38.8
VIS3648	36.0 in	48.0 in	38.8	50.8
VIS3660	36.0 in	60.0 in	38.8	62.8
M-VIS3636	900 mm	900 mm	38.8	38.8
M-VIS3648	900 mm	1200 mm	38.8	50.8
M-VIS3660	900 mm	1500 mm	38.8	62.8

