




OPTIC LOCKING SCREW 1/16 (1.5) ALLEN

2X ACTUATOR LOCKING SCREW 1/16 (1.5) ALLEN




	A
MAX	.27 (6.9)
MIN	.13 (3.3)



ACCEPTS .50 (12.7) DIA OPTICS


2X LOCATION HOLE, 5/64 DIA (2.00) PIN



2X 8-32 (M4) THD

.33 (8.4) CLEARS APERTURE AT 45°

45° ROTATED VIEW



9775(M) PINT-SIZE TOP ACTUATED MOUNT

 New Focus™