





# New Focus™

A Newport Corporation Brand

SCALE 1:2 [REV A]

8823-AC-UHV. MOUNT. 3 IN. AC. UHV. PICOMOTOR

©2009 NEWPORT CORP.. ALL RIGHTS RESERVED.



OPTIC RETAINER  
SCREW 5/64 HEX

4X PICOMOTOR  
KAPTON WIRE  
LENGTH 6.0 [152.4]

6X COUNTERBORE  
FOR 8-32 (M4)  
SOC HD CAP

NEGATIVE (-)  
POLARITY MARK

