



Newport®
Experience | Solutions

1830-R Series Quick Start Guide



Package Contents:

This package contains the following components:

- Model 1830-R or 1830-R-GPIB Power Meter
- Power Cord
- USB Flash Drive (contains Device User's Manual and Software)
- Calibration Certificate
- This Quick Start Guide

This package does not contain:

- Detector Devices, cabling or other accessories

Getting Started:

Before operating this device, refer to the Device User's Manual starting with the section on Safety Precautions. Verify the voltage select switch is set for your region per the User Manual.

The User Manual is contained on the USB flash drive enclosed in this package.

Software Information:

The root folder of the USB flash drive contains a readme file and the setup program to install the software. To install the USB driver and software double-click on **SETUP.EXE** file.



Verify voltage select switch is properly set for your region

Newport Corporation U.S.A.
Office: 800-222-6440

P/N: 90044528 REV. A