

2X THRU SLOTS FOR
1/4 (M6) MACHINE SCREWS

12X 8-32 (M4) THD

4X ALIGNMENT HOLES FOR
5/64 (2.0) ALLEN WRENCH

2.50 (63.5)

1.25 (31.8)

13X 1/4-20 (M6) THD

NOTE:
ENGLISH HOLES ON .50" SPACING
METRIC HOLES ON 12.5MM SPACING

4X 1/4-20 (M6) FOR
MOUNTING TRACK AT 45 DEG

1.75 (44.5)

3.50 (88.9)

.25 (6.4)

9021(M) BASE PLATE

 New Focus™