

.39 [9.9]
OPTIC
THICK

RETAINING RING
.035X40 THREAD
Ø .90 [22.9] CA

HOLDS Ø 1.00 [25.4]
OPTICS

.91 [23.1] .75 [19.1]

Ø .24 [6.0]
GUIDE ROD HOLE

2.10 [53.3]
NOM

1.75 [44.5]
5/64 [2.0]
HEX KEY

.16 [4.1]

.30 [7.6]

SECTION A-A

1.61 [41.0]
.48 [12.2]
1.96 [49.9]
NOM

