

IS- ϕ 1" Port
SMA-F.O Adaptor

IS3.3-Sub Assy


IS- ϕ 1" Port
Reducer To ϕ 0.5"

IS- ϕ 1.5" Port
Reducer To ϕ 1"

Front View

123 REF


120 REF



IS- ϕ 1" Port
 ϕ 1" Port plug

1-40 UNS-2A,
Detector Adaptor

918D-IG-OD1R, Sensor



ϕ 12.7

1/4-20 BSW ∇ 6 Min
Thread For Rod