



MODEL NUMBER	DESCRIPTION	LEVELING ARM	NO. OF BARB FITTING	FITTING CONNECTION
CM-225-1V	ISOLATOR, STANDARD	YES	1	TO AIR SUPPLY MANIFOLD
CM-225-1S	ISOLATOR, SECONDARY	NO	1	TO ISOLATOR, VALVED
	ISOLATOR, PRIMARY	YES	2	ONE TO ISOLATOR, SECONDARY ONE TO AIR SUPPLY MANIFOLD

MODEL NUMBER	DESCRIPTION
CM-225	SET OF FOUR ISOLATORS; TWO STANDARD, ONE PRIMARY, ONE SECONDARY
CM-225-3	SET OF THREE STANDARD ISOLATORS
UCM-225	SET OF FOUR CLEANROOM ISOLATORS: TWO STANDARD, ONE PRIMARY, ONE SECONDARY

- NOTES:
1. A STANDARD UNIT HAS ITS OWN VALVE TO CONTROL FLOATING.
  2. A SECONDARY UNIT IS PNEUMATICALLY HOOKED UP TO A PRIMARY UNIT, WHICH CONTROLS FLOATING FOR BOTH UNITS.
  3. DIMS IN [ ] ARE IN mm

